

***Biomarker komponensek szelektív kvantitatív
elemzése lovak szérum mintáinak
¹H-NMR spektrumában
a CRAFT szoftver segítségével***

Pelczer István

*Department of Chemistry
Frick Chemistry Laboratory
NMR Facility, Rm. B09
(ipelczer@princeton.edu)*

2

**Magyar NMR Munkabizottság
ülése**

Balatonszemes/Princeton
2025. május 8-9.




Esztergom

(my hometown in Hungary...)



Made by

Mrs. István Pelczer Sr.




*Hubert Sattler: Esztergom a Dunáról nézve, 1873
(View of Esztergom from the Danube)
Keresztény Múzeum, Esztergom, HU*

Akiknem az NMR-t köszönhetem...



The NMR Facility in the Frick Chemistry Laboratory



Co-workers



*Ken Conover
(NMR Specialist)*




*Sarah Ralston
Professor
Emerita (Rutgers U.)
and Fulbright
Scholar*



*Dr. Merrill Simpson '13
GH Cook Scholar
Rutgers University
Sechelt Animal Hospital, BC*

*The five + 1 adventurers... - Laboratory Learning
Program 2024 summer program, Princeton University*



Matheus Guedin

Christina Metaxas

James Chen

Ken Conover

Ayelet Yablon

Joseph Hu

Krish

Krishnamurthy

*ChemPacker LLC, San Jose, CA 95135
{CRAFT}*

Introduction

We have applied the CRAFT (Complete Reduction to Amplitude-Frequency Table) software tool to extract essential inflammatory biomarker signals quantitatively from the full 1D ¹H spectrum of horse serum samples. This simple approach is a competitive alternative to a variety of sophisticated experimental methods which aim to suppress the rest of the spectrum and keep only the peaks of interest [DIRE, JEDI]. The CRAFT selective extraction uses the full information content of the spectrum without any modifications or loss of sensitivity.

[DIRE] Lodge S, et al; Anal. Chem., 2021, 93, 3976–3986

[JEDI] Nitschke P, et al.; Anal. Chem., 2022, 94, 1333–1341



Chen J, Yablon A, Metaxas C, Guedin M, Hu J, Conover K, Simpson M, Ralston SL, Krishnamurthy K, Pelczer I;

“Suppress or not to suppress? – CRAFT IT! A targeted metabolomics case study extracting essential biomarker signals directly from the ^1H -NMR spectra of horse serum samples”

Metabolites (*accepted, in press*)


Osteochondrosis Dessicans, OCD

A growth-related defect, caused by abnormal calcification of cartilage matrices at the growth plates and articular surfaces.
Heritable, correlation between abnormal glucose/insulin metabolism and the presence of lesions reported.



*Also:
Obesity, Laminitis*

There are a few very specific inflammation markers in blood, namely lipoprotein N⁺Me signals of protein clusters (GlycA and GlycB)



Fuertes-Martín R, Correig X, Vallvé JC, Amigó N; Life 2021, 11, 1407

In addition, there is a composite resonance of phospholipids (SPC, 3.21 ppm), which is also an essential inflammatory biomarker. The relative integrals of these resonances provide clear indication of inflammatory conditions, often related to serious diseases like cancer or COVID-19 infection.


Supporting Information

J-Edited Diffusional (JEDI) Proton Nuclear Magnetic Resonance Spectroscopic Measurement of Glycoprotein and Supramolecular Phospholipid Biomarkers of Inflammation in Human Serum

Philipp Nitschke¹, Samantha Lodge¹, Torben Kimhofer¹, Reika Masuda¹, Sze-How Bong¹, Drew Hall¹, Hartmut Schäfer², Manfred Spraul², Nils Pompe², Tammo Diercks³, Ganeko Bernardo-Seisdedos³, José M Mato³, Oscar Millet³, Daniella Susic^{4,5}, Amanda Henry^{4,5}, Emad M El-Omar⁶, Elaine Holmes^{1,7}, John C. Lindon⁷, Jeremy K. Nicholson^{1,8}, and Julien Wist^{1,9}

Anal. Chem. 2022,
94, 1333–1341

Spectral comparison for 1D ¹H, CPMG, DIRE and JEDI



A : full spectrum

B : CPMG diffusion filtered

C : DIRE filtered


D : JEDI-PGPE

E : JEDI-PGSE

F : JEDI-PGDE

A representative $^1\text{H-NMR}$ spectrum of horse blood serum


(only water suppression applied)



Experimental

For this case study we used a cohort of 20 + 20 horse samples, healthy and OCD diagnosed, respectively, as described and discussed in the Cook Thesis of Merrill Simpson [George H. Cook Scholars Program Thesis, 2013, School of Environmental and Biological Sciences, Rutgers University]. The CRAFT analysis can be done retrospectively on many datasets acquired earlier as the water suppressed full spectrum is naturally part of the acquisition package.

The routine protocol – full spectrum with water suppression only



And now...

Krish Krishnamurthy


the magician of

CRAFT !

A schematic representation of CRAFT analysis workflow

(from MRC 2013, 51, 821–829)

CRAFT FID => Table => Spectrum




Extract information

- No precondition
- Iterative process
- Stop: Residual is random (i.e., no probability of exponential decay model)

Offset	Amplitude	decay rate	Phase
-3381.84	10.395	0.22727	95.0
-3071.89	200.091	0.21883	95.2
-2816.81	205.187	0.18838	90.9
-2685.30	44.507	0.10155	93.8
-2630.59	17.423	0.08195	96.5
-2498.65	110.320	0.13337	96.2
-2245.51	695.838	0.18692	94.7
.....			
-1698.98	20.776	0.14821	95.0
.....			
-1292.25	49.424	0.13773	93.5
-1162.14	46.537	0.13112	95.0
.....			
290.63	39.522	0.21040	89.4
2606.89	142.060	0.03783	90.1

Each row defined by a single exponential decay function



ppm	Amplitude	width	Phase
0.2286	10.395	1.40	95.0
0.8492	200.091	1.45	95.2
1.3600	205.187	1.69	90.9
1.6233	44.507	3.13	93.8
1.7329	17.423	3.88	96.5
1.9971	110.320	2.39	96.2
2.5040	695.838	1.70	94.7
.....			
3.5983	20.776	2.15	95.0
.....			
4.4127	49.424	2.31	93.5
4.6733	46.537	2.43	95.0
.....			
7.5823	39.522	1.51	89.4
12.2203	142.060	8.41	90.1


Interpretation

Sort
Combine
Filter
Search
Group
etc.


NMR rules
(User definable)

Each row defined by one or more exponential decay functions

A representative CRAFT model on a randomly selected spectrum from the set of 20+20 cohort




The SIMCA modelling: O-PLS-DA analysis of GlycB/SPC CRAFT amplitude ratios



Scores - CRAFT Data Tables - finalized data (2) (M4, OPLS-DA)

Colored according to classes in M4




Scaled proportionally to R2X, R2X[1] = 0.409, R2Xo[1] = 0.29,
Ellipse: Hotelling's T2 (95%)

SIMCA 18.0.1 - 2025/04/06 03:36:10 PM (UTC-4)

GlycA/GlycB

Scores - spcglycA (M3, OPLS-DA)

Colored according to classes in M3




Scaled proportionally to R2X, R2X[1] = 0.406, R2Xo[1] = 0.356,
Ellipse: Hotelling's T2 (95%)

SIMCA 18.0.1 - 8/6/2024 10:41:33 AM (UTC-4)

SPC/GlycA

Scores - CRAFT Data Tables - finalized data2 (M4, OPLS-DA)

Colored according to classes in M4



Scaled proportionally to R2X, R2X[1] = 0.546, R2Xo[1] = 0.454,
Ellipse: Hotelling's T2 (95%)

SIMCA 18.0.1 - 8/6/2024 10:36:36 AM (UTC-4)

SPC/(GlycA+GlycB)

Summary, conclusions, future work

We proved that *CRAFT is a viable and highly simplified tool* for diagnostics from a single and simple spectrum. This work is focusing on the three major inflammatory biomarker signals (GlycA, GlycB and SPC), but *the concept is valid to all identifiable biomarker signals* or patterns of those.


GlycA, GlycB and SPC can be quantitatively extracted from complex spectral environments, including overlap with other, irrelevant signals and broad background resonances.

This work proves that *OCD is an inflammatory condition indeed*, as the ratio of the GlycB and SPC amplitudes clearly separate the healthy and diseased cohorts using multivariate analysis, as shown by the O-PLS-DA results with very good validation. Other combinations of the selected signal amplitudes show limited or no good clustering.


CRAFT can be applied *RETROSPECTIVELY* !

We are in the process to extend this application to many other datasets collected earlier, including those on laminitis, age related studies, comparing metabolic signature of mustangs to that of domestic horses, falcon aspergillosis, etc.

Thank you!...



Scores - FINALFINAL (2) (M3, OPLS-DA)
Colored according to classes in M3



Scaled proportionally to R2X, R2X[1] = 0.382, R2Xo[1] = 0.249,
Ellipse: Hotelling's T2 (95%)



Check out our new ssNMR/DNP!...